en Instructions for use/Technical description

Support arm

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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



Support arm

Legend

- 1 Support arm
- 2 Proximal part of the support arm
- 3 Fastening element
- 4 Silicone washer
- 5 Clamping handle
- 6 Distal part of the support arm
- 7 Quick-release coupling
- 8 Sleeve
- 9 Instrument shaft

1. About this document

Note
General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply to the M-TRAC® support arm FF168R.
► For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING
Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION
Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The M-TRAC® support arm is used in minimally invasive surgery to support e.g. optical systems, cameras, instruments and trocars.

Additionally, a fastening element is used at the proximal end of the support arm and an adapter (e.g. NeuroPilot®) at the distal end.

The maximum support capacity of the support arm is 40 N.

2.1.2 Indications

Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Contraindications

No known contraindications.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to use, check that the product is in good working order.
- Keep the instructions for use accessible for the user.

Note
The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Sterility

Non-sterile packaged products

The product is delivered in an unsterile condition.

- Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Assembly

⚠ WARNING
Injuries caused by unstable fixation of optics, cameras, instruments, trocars, support arms and fastening elements!
► Apply appropriate care when supporting and fixing optical systems, cameras, instruments or trocars.
► Tighten the clamping handle and check the fixation of the support arm.
► Do not overload the support arm.
► Ensure that the instrument shaft engages properly.
► With the instrument shaft engaged, turn the sleeve clockwise.
► Firmly tighten the fastening element at the operating table rail and the support arm.

- To separate the sterile region from the unsterile region, slide silicone washer 4 over the rod, starting from the bottom end.
- Fasten holding arm 1 using fastening element 3 FF151R*, FF280R* or RT090R* to the operating table rail, see Fig. C.
*Optional accessory, not part of delivery scope
- Check for correct fixation of support arm 1: Slide the support arm 1 vertically and horizontally.
- Insert instrument shaft 9 into the sleeve 8, see Fig. D.
Quick-release coupling 7 is locked automatically.
- Position instrument shaft 9 at angles of 60°, see Fig. E.
- To secure quick-release coupling 7 against inadvertent opening, rotate sleeve 8 clockwise, see Fig. D.

2.4 Application

⚠ WARNING

Risk of injury and/or malfunction!

- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Always carry out a function test prior to each use of the product.

⚠ WARNING

Injuries caused by unstable fixation of optics, cameras, instruments or trocars!

- Apply appropriate care when supporting and fixing optical systems, cameras, instruments or trocars.
- Tighten the clamping handle and check the fixation of the support arm.
- Do not overload the support arm.

⚠ CAUTION

Damage to the joints of the support arm!

- Loosen the clamping handle before changing the position of the support arm.

Note

The holding arm 1 is fitted with insulation in the distal part 6 to isolate the connected instrument shaft 9 from earth.

- Ensure that the insulation during HF application is not bypassed.
- Lift the distal part of the support arm 6 and loosen the clamping handle 5 by turning it counterclockwise b, see Fig. A.
- Position support arm 1 in the operating field.
- Tighten clamping handle 5 by turning it clockwise a.
The joints of support arm 1 are fixated simultaneously.
- Check the fixation of support arm 1.

2.5 Disassembly

- Remove any optical systems, cameras, instruments or trocars coupled to the support arm from quick-release coupling 7.
- Turn sleeve 8 of quick-release coupling 7 counterclockwise.
Sleeve 8 of quick-release coupling 7 is unlocked.
- Pull sleeve 8 of quick-release coupling 7 towards instrument shaft 9, see Fig. B.
Instrument shaft 9 is pushed out from quick-release coupling 7.
- Remove fastening element 3 and silicone washer 4 from the support arm.

3. Validated reprocessing procedure

3.1 General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.ak-i.org, link to "AKI-Brochures", "Red brochure".

3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- Disassemble the product prior to cleaning, see Disassembly.
- Tighten the clamping handle by turning it clockwise.

3.6 Cleaning/Disinfection

3.6.1 Product-specific safety information on the reprocessing method

Danger to the patient! The product must not be manually reprocessed.

- ▶ Reprocess the product only in a mechanical washer/disinfector.

Damage to the clamping mechanism of the support arm during processing!

- ▶ Do not immerse the support arm in any fluid.
- ▶ Tighten the clamping handle by turning it clockwise.

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
 - that are approved for plastic material and high-grade steel,
 - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.

For products with plasma layers (e.g. Noir instruments), the layer is attacked or removed when special cleaning procedures with the addition of oxidizing chemicals (e.g. hydrogen peroxide H₂O₂) are used.

- ▶ Do not use oxidizing chemicals for cleaning.

3.6.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Mechanical neutral or mild-alkaline cleaning and thermal disinfection	■ Mount jaws protection on the product. ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ To flush the product: Use a flushing nozzle or flushing sleeve. ■ Keep working ends open for cleaning. ■ Place the product on the tray with all product links and joints open.	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical neutral or mild alkaline cleaning and thermal disinfecting

3.7 Mechanical cleaning/disinfection

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.7.1 Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	Neutral: ■ Concentrate: – pH neutral – <5 % anionic surfactant ■ 0.5 %* working solution Mildly alkaline: ■ Concentrate: – pH = 9.5 – <5 % anionic surfactant ■ 0.5 % solution
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner neutral

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.
- ▶ Repeat the cleaning/disinfecting process if necessary.

3.8 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.8.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.8.2 Functional test

Support arm malfunction due to improper maintenance!

- ▶ **Do not lubricate with oil the joints of the support arm.**
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.9 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.
- ▶ Open the clamping handle of the support arm by turning it counterclockwise as far as it will go.

3.10 Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.11 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Maintenance and service

4.1 Maintenance

This product is maintenance-free.
In cases of damage, stiffness or insufficient clamping force, please contact your national B. Braun/Aesculap agency, see Technical service.

4.2 Technical service

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ **Do not modify the product.**
- ▶ **For service and repairs, please contact your national B. Braun/Aesculap agency.**

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1601
Fax: +49 7461 16-2887
E-Mail: ats@aesculap.de
Other service addresses can be obtained from the address indicated above.

4.3 Accessories / spare parts

Art. no.	Designation
FF280R	Flexible fastening element with ball joint
RT090R	Flexible fastening element with tooth rim
FF151R	Rigid fastening element
RT046P	Adapter for endoscope holder
RT060R	NeuroPilot®
RT079R	Adapter for fixation on endoscope body/eyepiece
RT099R	MINOP® TREND adapter for holding arm
RT081R	Adapter for attaching RT055P
RT068R	MINOP® InVent adapter for holding arm

5. Disposal

Risk of infection due to contaminated products!

- ▶ **Adhere to national regulations when disposing of or recycling the product, its components and its packaging.**

Note
The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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