RespiChamber[®]*

HOSPITAL

DISPOSABLE VALVED HOLDING CHAMBER

Replacement copies of this insert can be obtained by visiting www.tmmml.com or by calling 1-800-265-5494

HEALTHCARE PROVIDER INFORMATION

The **RespiChamber[®]** HOSPITAL Disposable Valved Holding Chamber (VHC) is designed for short-term use under the supervision of a healthcare professional. Patients should visit their retail pharmacy to obtain an **AeroChamber Plus[®] Flow-Vu[®]** Chamber when they are discharged.

INTENDED USE: This chamber is intended to be used along with a metered dose inhaler (inhaler) to deliver aerosol medication to the patient's lungs.

This product is intended for single patient use within a hospital/clinic environment only.

BEFORE FIRST USE: Read instructions supplied with the inhaler and chamber.

NOTES

- Chamber should be discarded after approximately 7 days of use.
- Environmental conditions, storage and frequency of use can affect product lifespan.
- This product contains no latex.
- When not in use, store in suitable container.
- If you have any questions about the performance or usability of this product, please contact your TMMML Sales Representative.

CAUTIONS

- Do not leave the chamber unattended with children.
- Ensure directions for use have been read prior to use and are kept available at all times.

Trudell Medical International[®]

725 Third Street, London,
Ontario, Canada N5V 5G4
+ 519-455-7060
customerservice@trudellmed.com
www.trudellmed.com

Distributed by:
Trudell Medical Marketing Limited
758 Third Street, London,
Ontario, Canada N5V 5J7







+ 1-800-265-5494
tmmml@tmmml.com
www.tmmml.com

REGISTERED
ISO 13485
QUALITY SYSTEM

RespiChamber[®] HOSPITAL Disposable Chamber complies with CSA Standard Z264.1-02., Spacers and Holding Chambers for use with Metered Dose Inhalers.

DIRECTIONS FOR USE

<p>1 </p>	<p>Before use, carefully examine the chamber. Remove contamination if found. Replace immediately if any defect is noticed.</p>
<p>2 </p>	<p>Remove cap from inhaler.</p>
<p>3 </p>	<p>Shake the inhaler immediately before use as per the instructions supplied with it.</p>
<p>4 </p>	<p>Insert the inhaler into the back of the chamber.</p>
<p>5 </p>	<p>MOUTHPIECE: Put mouthpiece into mouth and ensure a secure seal. MASK: Slide mask onto the mouthpiece of the chamber. Align mask to inhaler per picture. Apply mask to face as shown and ensure an effective seal.</p>
<p>6 </p>	<p>MOUTHPIECE: Press the inhaler at the beginning of a slow inhalation. The patient should inhale slowly and deeply through the mouth until a full breath has been taken followed by a 5–10 second breath hold, if possible. Alternatively, the patient can keep their mouth tight on the mouthpiece and breathe slowly 2–3 times through the chamber after the inhaler is pressed. MASK: Press the inhaler at the beginning of a slow inhalation. Maintain seal for 5–6 breaths after inhaler is pressed.</p>
<p>7 </p>	<p>SLOW DOWN inhalation if the whistle sounds. It means the patient is inhaling too quickly.</p>
<p>8 </p>	<p>Follow instructions supplied with inhaler on the amount of time to wait before repeating steps 3–7. Administer one puff at a time.</p>

* Optional mask shown.

Product Number	Description	Case Size
T0110152800020	RespiChamber[®] HOSPITAL Disposable VHC	20
T0111353691050	Small Silicone Mask	50
T0111353791050	Medium Silicone Mask	50
T0111355791050	Large Silicone Mask	50

For home use prescribe **AeroChamber Plus[®] Flow-Vu[®]** Chamber for your metered dose inhaler patients.

- **AeroChamber Plus[®] Flow-Vu[®]** Chamber with Small Mask
0–18 months

- **AeroChamber Plus[®] Flow-Vu[®]** Chamber with Medium Mask
1–5 years

- **AC[®] Girlz[®]** Chamber with Mouthpiece
5 years +

- **AC[®] Boyz[®]** Chamber with Mouthpiece
5 years +

- **AeroChamber Plus[®] Flow-Vu[®]** Chamber with Mouthpiece

- **AeroChamber Plus[®] Flow-Vu[®]** Chamber with Large Mask

AeroChamber[®] VHC is covered by most insurance plans. A prescription is often required to obtain reimbursement.

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Trudell Medical International*

725 Third Street, London,
Ontario, Canada N5V 5G4

+ 519-455-7060

customerservice@trudellmed.com

www.trudellmed.com

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Trudell Medical Marketing Limited

758 Third Street, London,
Ontario, Canada N5V 5J7

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