**TRU TEST TISSUE ADHESIVE FLEXIBLE**  
**Package Insert**

Before using product, read the following information thoroughly.

**IMPORTANT!** This insert is designed to assist in using TruTest Tissue Adhesive Flexible. It is not a reference to surgical techniques.

**DESCRIPTION:**

TruTest Tissue Adhesive Flexible is a sterile, liquid topical skin adhesive composed of n-Butyl-2-Cyanoacrylate monomer with Myristic acid and Triacetin added to allow the film to remain flexible. TruTest Tissue Adhesive Flexible is colored with D&C Violet #2 dye in order to make it easier to see the thickness of the layer being applied. TruTest Tissue Adhesive Flexible is supplied in 0.5 ml single patient use plastic ampoules. Each ampoule is sealed within an aluminum pouch so the ampoule can remain sterile. TruTest Tissue Adhesive Flexible remains liquid until exposed to acidic, basic, alcohol, water or water-containing substances, including tissues. Tissue Adhesive Flexible cures (polymerizes exothermally) and forms a film that bonds to the underlying surface.

TruTest Tissue Adhesive Flexible may be used with an optional applicator tip. The applicator tip attaches to the open cannula of the TruTest Tissue Adhesive Flexible ampoule. The adhesive can then be applied to the skin surface.

*In vitro* studies have shown that TruTest Tissue Adhesive Flexible as a barrier to microbial penetration as long as the adhesive film remains intact. Clinical studies were not conducted to demonstrate microbial barrier properties and a correlation between microbial barrier properties and a reduction in infection have not been established. See *INSTRUCTIONS FOR USE* section for further information.

**INDICATIONS:**

TruTest Tissue Adhesive Flexible is intended for topical application to hold closed easily approximated skin edges of minimum-tension wounds from clean surgical incisions and simple, thoroughly cleansed, trauma-induced lacerations. TruTest Tissue Adhesive Flexible may be used in conjunction with, but not in place of, deep dermal sutures.

**CONTRAINDICATIONS:**

- TruTest Tissue Adhesive Flexible is not to be applied below the surface of the skin. The liquid adhesive will react exothermally with tissue; the polymerized adhesive is not absorbed by any tissues and may elicit a foreign body reaction.
- TruTest Tissue Adhesive Flexible is not to be applied to any internal organs, blood vessels, nerve tissue, mucosal surfaces or mucocutaneous junctions, areas with dense natural hair, or within the conjunctival sac of the eye.
- TruTest Tissue Adhesive Flexible is not to be applied to the surface of the eye. If the eyelids are accidentally bonded closed, release eyelashes with warm water by covering with a wet pad. The adhesive will bond to eye protein and will cause periods of weeping which will help to debond the adhesive. Keep the eye covered until debonding is complete – usually within 1 to 3 days. Do not force the eye open.

- TruTest Tissue Adhesive Flexible is not to be applied to wounds subject to high skin tension or on areas of increased skin tension such as the elbows, knees, or knuckles. TruTest Tissue Adhesive Flexible is not to be used in areas of skin excision.
- TruTest Tissue Adhesive Flexible is not to be applied to wounds that show evidence of infection, gangrene or wounds of decubitus etiology.
- TruTest Tissue Adhesive Flexible is not to be used on patients with known preoperative systemic infections, uncontrolled diabetes, or diseases or conditions that are known to interfere with the wound healing process.
- TruTest Tissue Adhesive Flexible is not to be used on patients with a known hypersensitivity to cyanoacrylate, formaldehyde, or the dye D&C Violet #2.

**WARNINGS:**

- TruTest Tissue Adhesive Flexible should be used only on wounds that have been thoroughly cleaned, debrided and have easily apposed wound edges.
- TruTest Tissue Adhesive generates a small amount of heat during polymerization and should not be applied to tissues that may be affected by such heat.
- TruTest Tissue Adhesive Flexible should always be applied very sparingly, either as minute drops or as a very thin film along the edges of the wound. Heavy application may cause thermal damage to tissues, and delayed healing may result.
- TruTest Tissue Adhesive Flexible should not be applied to wet wounds. Excess moisture, such as water or alcohol, may accelerate polymerization, resulting in the generation of excess heat.
- Use of TruTest Tissue Adhesive Flexible may result in localized sensitization or irritation reactions.
- Application and/or migration of TruTest Tissue Adhesive Flexible below the surface of the skin will impair the healing process by forming a barrier between tissue edges.
- TruTest Tissue Adhesive Flexible migration (leaking) below the epidermal surface may result in “tattooing” of the underlying tissue.
- TruTest Tissue Adhesive Flexible will readily adhere to most substrates. Care should be taken to avoid unwanted contact with the adhesive during polymerization. Polymerized adhesive can be removed from metal instruments with acetone. Accidental bonding of materials other than tissues may be reversed by peeling apart the adhered surfaces with the aid of warm soapy water, petroleum gel, saline solution, or 5% solution of sodium bicarbonate.
- Accidental bonding of unwanted skin may occur. Do not pull apart skin. Instead, accidental bonding of unintended areas of skin of the body can be corrected with the use of acetone or by soaking in warm water until the skin may be separated.

**PRECAUTIONS:**

- Wounds should be kept dry following closure with TruTest Tissue Adhesive Flexible. Do not apply topical medications following closure.
- In the event of spillage, TruTest Tissue Adhesive Flexible can be absorbed with talc. Flush area with water to solidify the adhesive.
- TruTest Tissue Adhesive Flexible has not been evaluated in patients with a history of hypertrophic scarring or keloid formation.

- Small quantities of TruTest Tissue Adhesive Flexible should be used during wound repair because use of excess TruTest Tissue Adhesive Flexible can result in tissue damage due to the cumulative heat dissipated during device polymerization.

**POSSIBLE ADVERSE EVENTS:**

Clinical experience with TruTest Tissue Adhesive Flexible suggests that the following adverse events may occur:

- Dehiscence with the need for re-closure
- Suspected infection
- Bonding to unintended tissues
- Thermal discomfort during polymerization
- Allergic reaction
- Foreign body reaction
- Tattooing
- Chronic non-healing of a wound

**HOW SUPPLIED:**

TruTest Tissue Adhesive Flexible is supplied in 0.5ml single patient use ampoules. Each ampoule is sealed within an aluminum pouch so the exterior of the ampoule can remain sterile. TruTest Tissue Adhesive Flexible is supplied in boxes of 10 individual ampoules.

**INSTRUCTIONS FOR USE:**

1. Inspect and clean the wound, provide local anesthesia for adequate cleansing and debridement of any devitalized structures, assure hemostasis, close the dermis as needed, and assure that surface edges are easily apposable before applying TruTest Tissue Adhesive Flexible.
2. For wounds at risk for tension, provide relief of potential stress along the wound line by approximating wound edges with subcuticular sutures before applying TruTest Tissue Adhesive Flexible to the skin surface.
3. Open the aluminum pouch to expose the sterile, single patient use ampoule.
4. The ampoule is held at its tip and rigorously shaken downwards to ensure that no adhesive remains within the cannula. Then the ampoule is opened by twisting off the ribbed tip of the cannula. It is advisable to hold the yellow part of the ampoule between two fingers, holding it vertically with the thin end pointing upward, while opening the ampoule. This will prevent any lost TruTest Tissue Adhesive Flexible from escaping from the ampoule when opening.
5. To express TruTest Tissue Adhesive Flexible from the ampoule tip, apply light pressure to the ampoule.
6. Appose the tissue edges and hold in apposition while applying TruTest Tissue Adhesive Flexible. Hold in apposition for approximately 60 seconds for TruTest Tissue Adhesive Flexible to polymerize and to prevent migration between wound edge.
7. Apply TruTest Tissue Adhesive Flexible to the easily apposed wound edges very sparingly, either as minute drops or as a very thin film along the top edges of the wound. Avoid heavy application.
8. After applying TruTest Tissue Adhesive Flexible, maintain light pressure along the wound line to maintain apposition for approximately 60 seconds to allow the adhesive to cure.
9. After TruTest Tissue Adhesive Flexible application is completed; discard the ampoule with any remaining adhesive by putting the ribbed tip over the thin end of the cannula.

**PATIENT INSTRUCTIONS:**

The following information should be shared with the patient:

- Patients may shower and/or bathe but allow only transient wetting of the treatment site. The site should not be soaked or exposed to prolonged wetness for 7-10 days or until polymerized film has sloughed off.
- Do not apply any medications or cream to the wound.
- Do not pull or pick at the wound.
- Avoid extreme physical activity that might dislodge or impact the wound surface.
- Report any unusual discomfort, redness, drainage, swelling and/or other concerns regarding your wound to your doctor.
- If they are recommended by your doctor, you should keep the wound dry with a non-stick, non-medicated and water-resistant bandage. While not required specifically for this product, your doctor in his/her expert opinion may recommend the use of a bandage as a special precaution for your individual circumstance.

**STORAGE:**

- TruTest Tissue Adhesive Flexible must be stored in its original sealed aluminum pouch.
- TruTest Tissue Adhesive Flexible should be kept in the aluminum pouch to maintain stability of the liquid adhesive at room temperature.
- TruTest Tissue Adhesive Flexible should be stored at ambient conditions in temperatures at or below 77°F (25° C) away from moisture and direct heat.
- Do not expose TruTest Tissue Adhesive Flexible to elevated temperatures (i.e., 104-140° F, 40-60° C) for more than eight hours.
- Do not use TruTest Tissue Adhesive Flexible after the expiration date shown on the aluminum pouch, preceded by the expiration symbol.

**STERILITY:**

TruTest Tissue Adhesive Flexible is packaged for single patient use. Do not re-sterilize. Do not use if package is opened or damaged.

**CUSTOMER SERVICE:**

For further information regarding TruTest Tissue Adhesive Flexible please contact Trudell Healthcare Solutions, Inc. at **1-800-265-5494** or **customerservice@trudellhs.com**.

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