

Reusable Scissors

Sharpness guaranteed




Aesculap Surgical Technologies - Endoscopy

AESCULAP[®]

*A member of the **B|BRAUN** group of companies.*








Reusable Scissors



Aesculap's Reposable Scissors System offers the convenience and guaranteed sharpness of a disposable scissor, with the cost savings of a reusable. The system provides:

- Superior handle-to-blade responsiveness
- Choices of jaw patterns and handle types
- Shaft insulated to the tip for patient and surgeon safety

Disposable Scissor Inserts – Reusable Quick-Snap™ Handle or Traditional Insulated Metal Handle

	Item No.	Description	Diameter	Length
	PO886	Disposable Hook Scissor Insert (10/box)	5 mm	31 cm
	PO887	Disposable Mini-Metz Scissor Insert (10/box)	5 mm	31 cm
	PO888	Disposable Metzenbaum Scissor Insert (10/box)	5 mm	31 cm
	PO889	Disposable Metzenbaum Scissor Insert (10/box)	5 mm	42 cm
	PO885SU	Disposable Slightly Curved Metzenbaum Scissor Insert (10/box)	5 mm	31 cm
	PO958R	Quick-Snap™ handle		
	PM950R	Traditional insulated metal handle		

Quick-Snap™ handle:

- Choice of durable PEEK™ or insulated metal
- Quick and easy assembly/disassembly

Assembly:




Expose red on rotator



Snap in disposable shaft

Disassembly:



Remove shaft by pulling back rotator

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2010 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com

DOC518 Rev. D 5M 7/10