

# For Dose Assurance Recommend *AeroChamber Plus<sup>®</sup> Flow-Vu<sup>®</sup>* aVHC

**NEW**

**Adult Small Mask**



Small Mask  
0-18 months

Medium Mask  
1-5 years

Mouthpiece

Adult  
Small Mask

Adult  
Large Mask

# Expertly Designed with Patients in Mind

## Adult Small Mask

Developed using anthropometric data analysis and refined through user and performance evaluations.

The new mask contours gently to the face with minimal effort, providing another facemask option for adult patients with smaller facial dimensions.

Facemasks are ideal for adult patients with co-morbidities or reduced cognitive ability or simply for increased ease of use.

Patients simply need to breathe tidally through the chamber to inhale the medication.

AeroChamber Plus* Flow-Vu* VHC	TMML Item #	Case Quantity	CMDL Number
Small Mask (orange)	T0110650300010	10	89428
Medium Mask (yellow)	T0110650200010	10	
Mouthpiece (blue)	T0110650500010	10	
<b>NEW</b> Adult Small Mask (purple)	T0110850800010	10	
Adult Large Mask (blue)	T0110650400010	10	



MD-171A-0618 - \* trade-marks and registered trade-marks of Trudell Medical International. † trade-marks and registered trade-marks of respective companies.  
© TMI 2018. All rights reserved. Manufactured in Canada with Canadian and USA Parts.



**TRUDELL**  
HEALTHCARE SOLUTIONS

758 Baransway Drive, London,  
Ontario, Canada, N5V 5G4

+1-800-265-5494

customerservice@trudellhs.com

www.trudellhs.com



**AeroChamber**  
Plus **Flow-Vu**

Anti-Static Valved Holding Chamber  
aerochamber.com