

**NEW**

# AeroChamber Plus Z STAT\*

Anti-Static Valved Holding Chamber



## Dose Assurance for Single Patient Use

The **AeroChamber Plus\* Z STAT\*** Anti-Static Valved Holding Chamber is a cost effective choice for hospitals who are looking for a high quality single patient use chamber.

The anti-static chamber provides consistent medication availability and can be used right out-of-package saving the hospital valuable time and resources.

A full range of products is available, designed specifically for each patient type.

Our brand of chambers is supported by numerous *in vitro* and *in vivo* studies and the safety and efficacy has been validated in third party clinical evaluations.



# Key features

## ComfortSeal\* Mask

(on mask models)

Contours gently to the face providing a **secure seal** and **comfortable fit**

Can be used with **all** commonly prescribed Metered Dose Inhalers



### EZ Flow Exhalation Valve

**Minimizes rebreathing** and directs exhaled flow away from the patient's face

### FLOWSignal Whistle (on adult models)

Alerts patients to **breathe slowly**

### Low resistance inhalation valve

Allows patients to **easily inhale** their medicine

### Anti-Static Chamber

Can be used directly out of package and provides **consistent medication** availability

The **AeroChamber Plus\* Z STAT\*** Anti-Static Valved Holding Chamber complies with CSA Standard Z264.1-02., Spacers and Holding Chambers for use with Metered Dose Inhalers.

	Reorder Number	Units per case
<b>AeroChamber Plus* Z STAT*</b> aVHC, Small Mask	T0112050300010	10
Medium Mask	T0112050200010	10
Large Mask	T0112050400010	10
Mouthpiece	T0112050100010	10



Designed and Developed in Canada

**HealthPRO**

Choice. Support. Results.

Please reference **04218** when placing your order.



**TRUDELL** HEALTHCARE SOLUTIONS

758 Baransway Drive, London, Ontario, Canada, N5V 5G4

+1-800-265-5494

customerservice@trudellhs.com

www.trudellhs.com



Manufactured by:

Trudell Medical International\*

725 Baransway Drive, London, Ontario, Canada, N5V 5G4

+1-866-510-0004

customerservice@trudellmed.com

www.trudellmed.com

ISO CLASS 8  
CLEANROOM MANUFACTURING