



Waterloo
Healthcare Medical Carts

PEDIATRIC CART



9-Drawer Cart, Color-coded for Emergency Pediatric Care

Our centralized lever-lock provides easy access for your tamper evident breakaway seals. It also serves as a "Transport Lock", keeping drawers in place during rapid transit.

DIMENSIONS

CART: 29"(w) x 24.5"(d) x 45"(h)
DRAWERS: 22"(w) x 16.5"(d)

Available in Elite Aluminum or Traditional Steel

ALUMINUM FEATURES

- Unicart Construction
- Self-closing Ball Bearing Slides
- Easy-roll Casters
- Dual Pullout Work Shelves
- Universal Accessory Mounting System
- Includes Drawer Divider Set
- Lightweight Construction

STEEL FEATURES

- Unicart Construction
- Self-closing Ball Bearing Slides
- Easy-roll Casters
- Single Pullout Work Shelf