

# Dual Cannula

The Dual Nasal Cannula is designed to provide oxygen to both nostrils and also sample CO<sub>2</sub> at the same time.

## Carbon Dioxide


CO<sub>2</sub> is sampled from both nostrils for side-stream CO<sub>2</sub> monitoring through the nares.

## Oxygen

Fresh oxygen is delivered to both nostrils via the ports.

## Soft Touch Nasal Cannula

Ears are the most common site for Medical Device Related Hospital Acquired Pressure Ulcers (MDR HAPU), with 35% of affected patients suffering tissue breakdown in this area<sup>1</sup>. Our nasal cannula range features Soft Touch tubing, minimizing friction and pressure on delicate tissues.



Male Luer



Female Luer



Microstream™ Luer

MALE LUER	FEMALE LUER	DESCRIPTION	TUBING	QTY
032-10-123U	032-10-124U	Adult Dual Nasal Cannula (Pigtail Luer)	7ft	25
032-10-125U	032-10-126U	Adult Dual Nasal Cannula	7ft	25
032-10-125AU	032-10-126AU	Adult Dual Nasal Cannula with Filter	7ft	25
032-10-127U	032-10-128U	Adult Dual Nasal Cannula	10ft	25
032-10-195U	032-10-193U	Adult Dual Oral Nasal Cannula	10ft	25
032-10-196U	032-10-194U	Adult Dual Oral Nasal Cannula with Filter	10ft	25



MICROSTREAM™ TYPE LUER



MICROSTREAM™ TYPE LUER + M/M ADAPTER

## Microstream™ Compatibility

Dual products with a Microstream™ type luer connector filter are compatible with capnography monitors and multiparameter monitors equipped with Microstream™ technology.

MICROSTREAM™ TYPE LUER	MICROSTREAM™ TYPE LUER + M/M ADAPTER	DESCRIPTION	TUBING	QTY
032-10-180U	032-10-180/1U	Adult Dual Oral Nasal Cannula with Filter	7ft	25
032-10-181U	032-10-181/1U	Adult Dual Oral Nasal Cannula with Filter	14ft	25
032-10-182U	032-10-182/1U	Adult Dual Nasal Cannula with Filter	7ft	25
032-10-183U	032-10-183/1U	Adult Dual Nasal Cannula with Filter	14ft	25
032-12-182U	032-12-182/1U	Pediatric Dual Nasal Cannula with Filter	7ft	25
032-12-183U	032-12-183/1U	Pediatric Dual Nasal Cannula with Filter	14ft	25
032-13-182U		Neonatal Dual Nasal Cannula with Filter	7ft	25



Universal Adapter

LUER	DESCRIPTION	TUBING	QTY
032-10-125CU	Adult Dual Nasal Cannula with Male Luer + Universal Adapter	7ft	25
032-12-125BU	Pediatric Dual Cannula with Male Luer + Slip Disc Filter	7ft	25