



Adult CO₂ monitoring mask

CE



Rx Only

**Flexicare Inc.**

15281 Barranca Parkway, Unit D,
Irvine, CA 92618
Toll Free: 1-800-985-3314
Tel: 949-450-9999
Fax: 949-450-9992
Email: inquiries@flexicare.com
www.flexicare.com

REF:





Flexicare Medical Ltd

Cynon Valley Business Park,
Mountain Ash. CF45 4ER. UK
Tel: +44 (0)1443 474647
Fax: +44 (0)1443 474222
Email: enquiries@flexicare.com
www.flexicare.com



5 055173 740159 >

YP5501-40U



Lpm/O ₂ % ≈					
(Lpm)	2	4	6	8	10
(O ₂ %)	28%	35%	40%	50%	60%

The British Thoracic Society recommends to use a minimum flow rate of 5 LPM with Medium Concentration Masks. However this depends on patient's inspiratory flow rate. All MC Masks are variable devices. It is recommended that for COPD patients a fixed concentration mask is used.