

# Endo-Bronch with Standard Tip (H-PSTDL)

## Double Lumen Endobronchial Tube, Left




### Clamp-less Y-Connector

- 90-degree valves eliminate the need for awkward clamps
- Pre-assembled for expediency
- Scope and suction catheter access in either Open (On) or Closed (Off) positions




### Low Friction Suction Catheters

- Textured outer wall for easier advancement
- Adjustable pressure suction connector



### Tear Resistant Cuff

- Durable PVC tracheal cuff resists tearing during insertion



### Anatomy Protection

- Low profile, smooth tracheal opening can reduce trauma during insertion
- More acute distal tip angle facilitates insertion into left main stem bronchus

# Ordering Information

Striving  
to be  
**#1 IN PATIENT  
SATISFACTION**

## Endo-Bronch with Standard Tip, Left-sided kit

### Kit includes:

- Left-sided Endo-Bronch tube with Standard Tip and pre-loaded stylet
- Pre-assembled, clamp-less Y-connector
- Two low friction suction catheters

**Not manufactured with  
natural rubber latex**

#### DESCRIPTION – 1 /PKG

#### PRODUCT NUMBER

Endo-Bronch with Standard Tip, Double Lumen Endobronchial Tube, Left, 28 Fr	H-PSTDL-28-1
Endo-Bronch with Standard Tip, Double Lumen Endobronchial Tube, Left, 32 Fr	H-PSTDL-32-1
Endo-Bronch with Standard Tip, Double Lumen Endobronchial Tube, Left, 35 Fr	H-PSTDL-35-1
Endo-Bronch with Standard Tip, Double Lumen Endobronchial Tube, Left, 37 Fr	H-PSTDL-37-1
Endo-Bronch with Standard Tip, Double Lumen Endobronchial Tube, Left, 39 Fr	H-PSTDL-39-1
Endo-Bronch with Standard Tip, Double Lumen Endobronchial Tube, Left, 41 Fr	H-PSTDL-41-1

**For information or product support**  
**1-800-421-0024 ■ [customercare@salterlabs.com](mailto:customercare@salterlabs.com)**

2365 Camino Vida Roble – Carlsbad, CA 92011  
Phone: 661-854-3166, 1-800-421-0024, 1-800-235-4203  
Fax: 661-854-3850  
Toll-Free Fax: 1-800-628-4690  
E-mail: [customercare@salterlabs.com](mailto:customercare@salterlabs.com)  
[www.salterlabs.com](http://www.salterlabs.com)

SLML-174 Rev A April 2017



**SALTER LABS®**  
DESIGN INNOVATIONS  
FOR RESPIRATORY CARE