

- Sterile
- Single-Use
- DEHP Free
- Clear Indicators



LarySeal Pro aids rapid and secure airway management, improves patient safety with gastric access to reduce risk of pulmonary aspiration, and facilitates endotracheal intubation during emergent situations.

Reduced Risk of Aspiration

An integrated suction catheter port for efficient removal of fluids and gastric contents.

Directed Intubation

The guide re-directs the tip of the endotracheal tube more anterior towards the trachea improving placement.

Reinforced Bite Guard

Minimizing risk of pressure to soft tissues and protects from occlusion if patient bites down while also reducing likelihood of dental damage.

Epiglottic Flap

Integral fenestrated flap protects from blockage with minimal increase in flow resistance, and elevates the epiglottis for easy, trauma-free ET tube or bronchoscope insertion.



PART NUMBER	DESCRIPTION	MAX ET TUBE SIZE	MAX SUCTION CATHETER SIZE	COLOR	PATIENT WEIGHT RANGE	QTY
038-94-810U	LarySeal Pro Size 1	3.5	8Fr	○	<5kg	10
038-94-815U	LarySeal Pro Size 1.5	4	10Fr	●	5kg - 10kg	10
038-94-820U	LarySeal Pro Size 2	5	12Fr	●	10kg - 20kg	10
038-94-825U	LarySeal Pro Size 2.5	5.5	14Fr	●	20kg - 30kg	10
038-94-830U	LarySeal Pro Size 3	7.5	16Fr	●	30kg - 50kg	10
038-94-840U	LarySeal Pro Size 4	8	18Fr	●	50kg - 70kg	10
038-94-850U	LarySeal Pro Size 5	8	18Fr	●	70kg - 100kg	10

LarySeal™

A comprehensive range of Laryngeal Mask Airways that adapt to the contour of the oropharyngeal airway for easy first time insertion success. Each of the LarySeal types share multiple common features and are manufactured to give you clinical confidence.

Crush Resistant Tube

The thickness of the main tube wall reduces risk of crushing, while the star lumen ensures maintained patency.

Secure Seal for Airway Management

The cuff of the mask forms itself to the contours of the oropharyngeal area, providing a secure seal.

Easy Indication of the Cuff Pressure through Pilot Balloon

Thin wall pilot balloon indicates the inflation of the cuff.

Universal Connector


15mm Male connector provides a universal connection.

Color-Coded Pilot Balloon

Provides a visual confirmation of the size as well as corresponding to the cuff inflator.

Clear View

Deep and open port reduces risk of blockage by epiglottis and unrestricted access.




LarySeal Blue

- Medical grade silicone
- Single-use & sterile
- Eco-friendly slimline pack option



LarySeal Clear

- Satin finished PVC cuff
- Single-use & sterile
- Clear Plus option includes lubricating jelly and the LarySeal Cuff Inflator



LarySeal Multiple


- Medical grade silicone
- Autoclavable up to 40 cycles
- Durable, high temperature resistant material



LarySeal Flexi

- Satin finished PVC cuff
- Flexible main tube
- Single-use & sterile
- Reinforced kink-free construction

LARYSEAL MULTIPLE	LARYSEAL BLUE	LARYSEAL CLEAR	LARYSEAL CLEAR PLUS	LARYSEAL FLEXI		COLOR	PATIENT WEIGHT RANGE	QTY	LARYSEAL BLUE SLIMLINE PACK	QTY
038-94-110U	038-94-210U	038-94-310U	038-94-310P		LarySeal Size 1	○	<5kg	10	038-94-610U	20
038-94-115U	038-94-215U	038-94-315U	038-94-315P		LarySeal Size 1.5	●	5kg - 10kg	10	038-94-615U	20
038-94-120U	038-94-220U	038-94-320U	038-94-320P	038-94-520U	LarySeal Size 2	●	10kg - 20kg	10	038-94-620U	20
038-94-125U	038-94-225U	038-94-325U	038-94-325P	038-94-525U	LarySeal Size 2.5	●	20kg - 30kg	10	038-94-625U	20
038-94-130U	038-94-230U	038-94-330U	038-94-330P	038-94-530U	LarySeal Size 3	●	30kg - 50kg	10	038-94-630U	20
038-94-140U	038-94-240U	038-94-340U	038-94-340P	038-94-540U	LarySeal Size 4	●	50kg - 70kg	10	038-94-640U	20
038-94-150U	038-94-250U	038-94-350U	038-94-350P	038-94-550U	LarySeal Size 5	●	70kg - 100kg	10	038-94-650U	20



LarySeal Cuff Inflator

- Color coded to match LarySeal's pilot balloon
- Quick identification of fill volume

		QTY
038-94-001U	LarySeal Color Coded Cuff Inflator	10