




The
right
connection
is **everything.**

Medical Fittings

Product Catalog



Medical Fittings Product Catalog

Overview.....	4
Ohmeda Style Quick Connect.....	6
Ohmeda Style Winged Quick Connect.....	7
Compact Ohmeda Style Check Unit (Female Coupler).....	8
Ohmeda Style Check Unit (Female Coupler).....	9
Chemetron Style Quick Connect.....	10
Chemetron Style Check Unit (Female Coupler)	11
DISS Female.....	12
DISS Male.....	14
Puritan-Bennett Style Quick Connect	16
Puritan-Bennett Style Check Unit (Female Coupler)	17
Schrader Style Quick Connect	18
Schrader Style Check Unit (Female Coupler).....	19
Medstar Style Quick Connect.....	20
Oxequip Style Quick Connect	21
“Y” Block Assemblies and DISS “Y” Assemblies	22
“T” Bar Assemblies and DISS “T” Assemblies	23
T-Branch, Y-Block, Hex/Quad Block and Rail Blocks.....	26
Extension Manifolds, Pole Clamp and Ferrules.....	27
Medical Grade Hose	28
Tubing Nipples	29
NPT Pipe Fittings	30
Condensation Traps.....	32
Neb Drives and Rescue Flows.....	33
International Fittings.....	34

The right connection is everything.

Quality products to make the right connection.

Precision Medical is a global pioneer in respiratory care focused on developing and manufacturing high-quality medical devices. With three main product lines spanning the hospital and homecare markets, our product portfolio covers infants through advanced age and ranges from high-quality oxygen administration to aspiration and suction products.

Precision Medical's durable medical fittings, which are made in the USA, allow full customization and targeting to meet the needs of our customers and users. Our fittings are found in hospitals, surgery centers, and long-term/sub-acute facilities worldwide. We manufacture adapters to fit every need such as both male and female couplers and quick connects for all major styles including Ohmeda, Chemetron, Puritan-Bennett, Schrader, and DISS connections. Customized configurations and hose assemblies are available to meet your unique requirements.

At Precision Medical, we take great pride in providing exceptional customer service. Our knowledgeable representatives are here to help you match the right product to your specific need and have it shipped quickly. We are 100% committed to providing products that are safe, accurate, durable and easy to maintain.

Precision Medical. The Difference is Precision.

Ohmeda Style Quick Connect

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Male	2101	2201	2301	2401	2601
 1/8" NPT Male x 1.5" Long			2301X		
 1/8" NPT Female	2102	2202	2302	2402	2602
 1/4" NPT Male	2103	2203	2303	2403	2603
 DISS Male (O2)	2105				
 DISS Male		2205	2305	2405	2605
 DISS Hex Nut	2106	2206	2306	2406	2606
 DISS Hand Tight	2108	2208	2308	2408	2608
 Hose Barb	1/4" Hose Barb 2111	1/4" Hose Barb 2211	1/4" Hose Barb 2311 5/16" Hose Barb 2312	1/4" Hose Barb 2411	1/4" Hose Barb 2611 5/16" Hose Barb 2615
 Hose Barb Swivel	1/4" Hose Barb 2111S	1/4" Hose Barb 2211S	1/4" Hose Barb 2311S 5/16" Hose Barb 2312S	1/4" Hose Barb 2411S	1/4" Hose Barb 2611S 5/16" Hose Barb 2615S
 DISS Male with Check Valve (O2)	2114				
 DISS Male with Check Valve		2214	2314	2414	2614


Ohmeda Style Winged Quick Connect

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 1/8" NPT Male	2101W	2201W	2301W	2401W
 1/8" NPT Male x 1.5" Long			2301WX	
 1/8" NPT Female	2102W	2202W	2302W	2402W
 1/4" NPT Male	2103W	2203W	2303W	2403W
 DISS Male (O2)	2105W			
 DISS Male		2205W	2305W	2405W
 DISS Hex Nut	2106W	2206W	2306W	2406W
 DISS Hand Tight	2108W	2208W	2308W	2408W
 Hose Barb	1/4" Hose Barb 2111W	1/4" Hose Barb 2211W	1/4" Hose Barb 2311W 5/16" Hose Barb 2312W	1/4" Hose Barb 2411W
 Hose Barb Swivel	1/4" Hose Barb 2111WS	1/4" Hose Barb 2211WS	1/4" Hose Barb 2311WS 5/16" Hose Barb 2312WS	1/4" Hose Barb 2411WS
 DISS Male with Check Valve (O2)	2114W			
 DISS Male with Check Valve		2214W	2314W	2414W

Compact Ohmeda Style Check Unit
(Female Coupler)

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Female	2160	2260	2360	2460	2660
 1/8" NPT Male	2161	2261	2361	2461	2661
 1/4" NPT Male	2163	2263	2363	2463	2663
 Hose Barb	1/4" Hose Barb 2164	1/4" Hose Barb 2264	1/4" Hose Barb 2364 5/16" Hose Barb 2367	1/4" Hose Barb 2464	1/4" Hose Barb 2664 5/16" Hose Barb 2667
 DISS Male (O2)	2165				
 DISS Male		2265	2365	2465	2665
 DISS Hex Nut	2166	2266	2366	2466	2666
 DISS Hand Tight	2168	2268	2368	2468	2668

Ohmeda Style Check Unit
(Female Coupler)

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Female	2150	2250	2350	2450	2650
 1/8" NPT Male	2151	2251	2351	2451	2651
 1/4" NPT Male	2153	2253	2353	2453	2653
 Hose Barb	1/4" Hose Barb 2154	1/4" Hose Barb 2254	1/4" Hose Barb 2354 5/16" Hose Barb 2357	1/4" Hose Barb 2454	1/4" Hose Barb 2654 5/16" Hose Barb 2657
 DISS Male (O2)	2155				
 DISS Male		2255	2355	2455	2655
 DISS Hex Nut	2156	2256	2356	2456	2656
 DISS Hand Tight	2158	2258	2358	2458	2658

Chemetron Style Quick Connect



	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Male	3101	3201	3301	3401	3601
 1/8" NPT Male x 1.5" Long			3301X		
 1/8" NPT Female	3102	3202	3302	3402	3602
 1/4" NPT Male	3103	3203	3303	3403	3603
 DISS Male (O2)	3105				
 DISS Male		3205	3305	3405	3605
 DISS Hex Nut	3106	3206	3306	3406	3606
 DISS Hand Tight	3108	3208	3308	3408	3608
 Hose Barb	1/4" Hose Barb 3111	1/4" Hose Barb 3211	1/4" Hose Barb 3311 5/16" Hose Barb 3312	1/4" Hose Barb 3411	1/4" Hose Barb 3611 5/16" Hose Barb 3615
 Hose Barb Swivel	1/4" Hose Barb 3111S	1/4" Hose Barb 3211S	1/4" Hose Barb 3311S 5/16" Hose Barb 3312S	1/4" Hose Barb 3411S	1/4" Hose Barb 3611S 5/16" Hose Barb 3615S
 DISS Male with Check Valve (O2)	3114				
 DISS Male with Check Valve		3214	3314	3414	3614



Chemetron Style Check Unit
(Female Coupler)

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Female	3150	3250	3350	3450	3650
 1/8" NPT Male	3151	3251	3351	3451	3651
 1/4" NPT Male	3153	3253	3353	3453	3653
 Hose Barb	1/4" Hose Barb 3154	1/4" Hose Barb 3254	1/4" Hose Barb 3354 5/16" Hose Barb 3357	1/4" Hose Barb 3454	1/4" Hose Barb 3654 5/16" Hose Barb 3657
 DISS Male (O2)	3155				
 DISS Male		3255	3355	3455	3655
 DISS Hex Nut	3156	3256	3356	3456	3656
 DISS Hand Tight	3158	3258	3358	3458	3658

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Dioxide
 Hex Nut- 1/8" NPT with 1" Stem	0106						
 Hex Nut-1/8" NPT with 1" Stem Swivel	0106S						
 Hex Nut- 1/8" NPT with 1½" Stem	0107	0207	0307	0407	0507		504248
 Hex Nut- 1/4" NPT with 1½" Stem					0506	0606	
 Hex Nut- Hose Barb	1/4" Hose Barb 0115	1/4" Hose Barb 0216	1/4" Hose Barb 0316 5/16" Hose Barb 0317	1/4" Hose Barb 0416	1/4" Hose Barb 0516	1/4" Hose Barb 0616 5/16" Hose Barb 0617	1/4" Hose Barb 504154
 Hex Nut- Hose Barb Swivel with 1½" Stem	1/4" Hose Barb 0115S	1/4" Hose Barb 0216S	1/4" Hose Barb 0316S 5/16" Hose Barb 0317S	1/4" Hose Barb 0416S	1/4" Hose Barb 0516S	1/4" Hose Barb 0616S 5/16" Hose Barb 0617S	1/4" Hose Barb 504154S
 Oxygen DISS Female x Oxygen DISS Female	0104						

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD	Carbon Dioxide
 Hand Tight- 1/8" NPT with 1" Stem	0110						
 Hand Tight- Hose Barb Swivel	0110S						
 Hand Tight- 1/8" NPT with 1½" Stem	0111	0209	0309	0409	0509		507883
 Hand Tight- 1/4" NPT with 1½" Stem					0508	0608	
 Hand Tight- Hose Barb	1/4" Hose Barb 0118	1/4" Hose Barb 0218	1/4" Hose Barb 0318 5/16" Hose Barb 0319	1/4" Hose Barb 0418	1/4" Hose Barb 0518	1/4" Hose Barb 0618 5/16" Hose Barb 0619	1/4" Hose Barb 504287
 Hand Tight- Hose Barb Swivel with 1½" Stem	1/4" Hose Barb 0118S	1/4" Hose Barb 0218S	1/4" Hose Barb 0318S 5/16" Hose Barb 0319S	1/4" Hose Barb 0418S	1/4" Hose Barb 0518S	1/4" Hose Barb 0618S 5/16" Hose Barb 0619S	1/4" Hose Barb 504287S

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 Nitrogen	 WAGD	 Carbon Oxide
 DISS Male to 1/8" NPT Male	0701	0731	0751	0771	0791		
 DISS Male to 1/4" NPT Male	0703	0733	0753	0773	0792	0794	8946
 1/4" Hose Barb	0709					0601 Demand Check 0602	
 1/4" Hose Barb with Knob	0710	0740	0760	0780	0795		
 1/4" Hose Barb with Knob & Check	0711	0741	0761	0781	0796		
 1/4" Hose Barb with Knob Swivel	0710S	0740S	0760S	0780S	0795S		
 1/4" Hose Barb with Knob & Check Swivel	0711S	0741S	0761S	0781S	0796S		


	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 Nitrogen	 WAGD	 Carbon Oxide
Demand Check Valve							
 DISS Male with Demand Check Valve to 1/8" NPT Male	0712	0742		0782	0797		
 DISS Male with Demand Check Valve to 1/4" NPT Male	0713	0743	0763	0783	0798	0799	8947

One Way Check Valve

	O2, DISS Male x 1/4" NPT Male	508373
	Air, DISS Male x 1/4" NPT Male	508374
	1/4" NPT Female x 1/4" NPT Male	508375

 Oxygen DISS 90° Elbow Male x Hex Nut	1" Stem 0705 1½" Stem 0706	 Oxygen Power Take-Off (1/8" NPT)	0150
 Oxygen DISS Male x DISS Male	0704	 Air Power Take-Off (1/8" NPT)	0151
 Oxygen DISS Male x 1/8" NPT Female	0702		

Puritan-Bennett Style Quick Connect

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD
 1/8" NPT Male	4101	4201	4301	4401	4601
 1/8" NPT Male x 1.5" Long			4301X		
 1/8" NPT Female	4102	4202	4302	4402	4602
 1/4" NPT Male	4103	4203	4303	4403	4603
 DISS Male with Knob (O2)	4105				
 DISS Male with Knob		4205	4305	4405	4605 (no knob)
 DISS Hex Nut	4106	4206	4306	4406	4606
 DISS Hand Tight	4108	4208	4308	4408	4608
 Hose Barb with Knob	1/4" Hose Barb 4111	1/4" Hose Barb 4211	1/4" Hose Barb 4311 5/16" Hose Barb 4312 (no knob)	1/4" Hose Barb 4411	1/4" Hose Barb 4611 (no knob) 5/16" Hose Barb 4615 (no knob)
 Hose Barb Swivel	1/4" Hose Barb 4111S	1/4" Hose Barb 4211S	1/4" Hose Barb 4311S 5/16" Hose Barb 4312S	1/4" Hose Barb 4411S	1/4" Hose Barb 4611S 5/16" Hose Barb 4615S
 DISS Male with Knob and Check Valve (O2)	4114				
 DISS Male with Knob and Check Valve		4214	4314	4414	4614 (no knob)

Puritan-Bennett Style Check Unit
(Female Coupler)

	 Oxygen	 Air	 Vacuum
 1/8" NPT Female	4150	4250	4350
 1/8" NPT Male	4151	4251	4351
 1/4" NPT Male	4153	4253	4353
 Hose Barb	1/4" Hose Barb 4154	1/4" Hose Barb 4254	1/4" Hose Barb 4354 5/16" Hose Barb 4357
 DISS Male (O2)	4155		
 DISS Male		4255	4355
 DISS Hex Nut	4156	4256	4356
 DISS Hand Tight	4158	4258	4358

Schrader Style Quick Connect

Schrader Style Check Unit
(Female Coupler)


	Oxygen	Air	Vacuum	Nitrous Oxide
	6101	6201	6301	6401
1/8" NPT Male				
			6301X	
1/8" NPT Male x 1.5" Long				
	6102	6202	6302	6402
1/8" NPT Female				
	6103	6203	6303	6403
1/4" NPT Male				
	6105			
DISS Male (O2)				
		6205	6305	6405
DISS Male				
	6106	6206	6306	6406
DISS Hex Nut				
	6108	6208	6308	6408
DISS Hand Tight				
	1/4" Hose Barb 6111	1/4" Hose Barb 6211	1/4" Hose Barb 6311	1/4" Hose Barb 6411
Hose Barb				
	6114			
DISS Male and Check Valve (O2)				
		6214	6314	6414
DISS Male and Check Valve				
Schrader Piston Identification				
	Hex Piston (standard)		Swivel Piston (add "S" after part number)	

	Oxygen	Air	Vacuum	Nitrous Oxide	Air/Nitrogen
	6150	6250	6350	6450	
1/8" NPT Female					
	6151	6251	6351	6451	
1/8" NPT Male					
	6153	6253	6353	6453	
1/4" NPT Male					
	1/4" Hose Barb 6154	1/4" Hose Barb 6254	1/4" Hose Barb 6354	1/4" Hose Barb 6454	
Hose Barb					
	6155				
DISS Male (O2)					
		6255	6355	6455	
DISS Male					
	6156	6256	6356	6456	Schrader Air Coupler x Nitrogen DISS 6556
DISS Hex Nut					
	6158	6258	6358	6458	Schrader Air Coupler x Nitrogen DISS 6558
DISS Hand Tight					

Medstar Style Quick Connect

	 Oxygen	 Air	 Vacuum
 1/8" NPT Male	8101	8201	8301
 1/8" NPT Male x 1.5" Long			8301X
 1/8" NPT Female	8102	8202	8302
 1/4" NPT Male	8103	8203	8303
 DISS Male (O2)	8105		
 DISS Male		8205	8305
 DISS Hex Nut	8106	8206	8306
 DISS Hand Tight	8108	8208	8308
 Hose Barb	1/4" Hose Barb 8111	1/4" Hose Barb 8211	1/4" Hose Barb 8311
 Hose Barb Swivel	8111S	8211S	8311S
 DISS Male with Check Valve (O2)	8114MS		
 DISS Male with Check Valve		8214	8314

Oxequip Style Quick Connect

	 Oxygen	 Air	 Vacuum
 1/8" NPT Male	5101	5201	5301
 1/8" NPT Male x 1.5" Long			5301X
 1/8" NPT Female	5102	5202	5302
 1/4" NPT Male	5103	5203	5303
 DISS Male with Knob	5105		
 DISS Male with Knob		5205	5305
 DISS Hex Nut	5106	5206	5306
 DISS Hand Tight	5108	5208	5308
 Hose Barb	1/4" Hose Barb 5111	1/4" Hose Barb 5211	1/4" Hose Barb 5311
 DISS Male with Knob and Check Valve	5114		
 DISS Male with Knob and Check Valve		5214	5314

“Y” Block Assemblies and DISS “Y” Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen
 Ohmeda Style	2152Y	2252Y	2352Y	2452Y	
with Winged Quick Connect	2152WY	2252WY	2352WY	2352WY	
 Compact Ohmeda Style	2162Y	2262Y	2362Y	2462Y	
with Winged Quick Connect	2162WY	2262WY	2362WY	2462WT	
 Chemetron Style	3152Y	3252Y	3352Y	3452Y	
 Puritan-Bennett Style	4152Y	4252Y	4352Y		
 Schrader Style	6152Y	6252Y	6352Y	6452Y	
 DISS “Y” Assemblies	7211	7212	7213	7214	7240

“T” Bar Assemblies and DISS “T” Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide
 Ohmeda Style	2152T	2252T	2352T	2452T
with Winged Quick Connect	2152WT	2252WT	2352WT	2452WT
 Ohmeda Style Offset Right	2152T-R	2252T-R	2352T-R	2452T-R
with Winged Quick Connect	2152WT-R	2252WT-R	2352WT-R	2452WT-R
 Ohmeda Style Offset Left	2152T-L	2252T-L	2352T-L	2452T-L
with Winged Quick Connect	2152WT-L	2252WT-L	2352WT-L	2452WT-L
 Compact Ohmeda Style	2162T	2262T	2362T	2462T
with Winged Quick Connect	2162WT	2262WT	2362WT	2462WT
 Compact Ohmeda Style Offset Right	2162T-R	2262T-R	2362T-R	2462T-R
with Winged Quick Connect	2162WT-R	2262WT-R	2362WT-R	2462WT-R
 Compact Ohmeda Style Offset Left	2162T-L	2262T-L	2362T-L	2462T-L
with Winged Quick Connect	2162WT-L	2262WT-L	2362WT-L	2462WT-L

“T” Bar Assemblies and DISS “T” Assemblies


	Oxygen	Air	Vacuum	Nitrous Oxide
 Chemetron Style	3152T	3252T	3352T	3452T
 Chemetron Style Offset Right	3152T-R	3252T-R	3352T-R	3452T-R
 Chemetron Style Offset Left	3152T-L	3252T-L	3352T-L	3452T-L
 Puritan-Bennett Style	4152T	4252T	4352T	
 Puritan-Bennett Style Offset Right	4152T-R	4252T-R	4352T-R	
 Puritan-Bennett Style Offset Left	4152T-L	4252T-L	4352T-L	

“T” Bar Assemblies and DISS “T” Assemblies

	Oxygen	Air	Vacuum	Nitrous Oxide
 Schrader Style	6152T	6252T	6352T	6452T
 Schrader Style Offset Right	6152T-R	6252T-R	6352T-R	6452T-R
 Schrader Style Offset Left	6152T-L	6252T-L	6352T-L	6452T-L
 DISS “T” Assemblies	7111T	7112T		7114T
 DISS “T” Assemblies Offset Right	7111T-R	7112T-R		7114T-R
 DISS “T” Assemblies Offset Left	7111T-L	7112T-L		7114T-L

T-Branch, Y-Block, Hex/Quad Block and Rail Blocks

			
T-Branch	Y-Block	Hex Block (3 outlet)	Quad Block (4 outlet)
1/8" NPT Standard 7100	1/8" NPT Female 7201	1/8" NPT Female 9054	
1/8" NPT Offset 7100L/R	1/4" NPT Female 7200	1/4" NPT Female 9052	1/4" NPT Female 9051

	Oxygen	Air	Vacuum
	7217	7218	7219
Hex Block (1/8" NPT Female Inlet)			
	7220	7221	7222
Quad Block Assembly (1/8" NPT Female Inlet)			
Rail Block Assemblies All Rail Block Assemblies include Wall Bracket and Plate	7223	7224	7225
 Plate Wall Bracket Rail Block Assembly (DISS Male Inlet)	Replacement Parts 7226 — Mounting Bracket 9056 — Mounting Plate 9054 — Hex Block, 3 outlets 9058 — Screw (2 required)		

Custom assemblies available upon request.

Extension Manifolds, Pole Clamp and Ferrules

Extension Manifold

FEATURES:
Economical expansion of single wall outlet.
Permanently mounts to wall.
2 foot conductive hose with 1/8" NPT Female inlet is standard.
Outlets are DISS Male with Check Valves.

			
	Oxygen	Air	Vacuum
	7230	7233	
Extension Manifold-2 Outlet			
			7236
Extension Manifold-2 Outlet with Swivel			
	7231	7234	
Extension Manifold-3 Outlet			
	7232	7235	
Extension Manifold-4 Outlet			
	505568		
Oxygen Manifold			

Pole Clamp

FEATURES:
1/8" NPT Female holes (top and front) - Comes with 1 plug on bottom

	7227A shown	1" Pole 7216
	7227B shown	1-1/2" Pole 7227
Pole Clamp		

Ferrule

	1/4" Hose		5/16" Hose	
	3/4" Long	1" Long	3/4" Long	1" Long
	9030	507712	503780	507501
Ferrule	Add -25 or -50 at end of Part Number for those quantities			

FEATURES:

Kink resistant and light weight.

FDA medical grade conductive hose.

Reinforced with tough nylon braiding.

Conveniently packaged in 250 foot rolls or as noted.

Color coded to conform to CGA-C9 standards.

Maximum working pressure to 200 psi @ 70°F.

Meets NFPA 99 requirements for conductivity (conductive hose only).

Keep out of sunlight and in a climate controlled atmosphere

Non-Conductive Hose			
Oxygen	Green		(1/4" ID-250' Roll) 9020-250
Air	Yellow		(1/4" ID-250' Roll) 9021-250
Vacuum	White		(1/4" ID-250' Roll) 9022-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll) 9023-250
Nitrogen	Black		(1/4" ID-250' Roll) 9028-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll) 9033-250
WAGD	Purple		(5/16" ID-250' Roll) 9038-250
	Clear		(1/4" ID-100' Roll) 504202-100

Conductive Hose			
Oxygen	Green		(1/4" ID-250' Roll) 9120-250
Air	Yellow		(1/4" ID-250' Roll) 9121-250
Vacuum	White		(1/4" ID-250' Roll) (5/16" ID-250' Roll) 9122-250 9135-250
Nitrous Oxide	Blue		(1/4" ID-250' Roll) 9123-250
Nitrogen	Black		(1/4" ID-250' Roll) 9128-250
Carbon Dioxide	Gray		(1/4" ID-250' Roll) 9125-250
WAGD	Purple		(1/4" ID-250' Roll) 9124-250
ISO Air	Black & White Striped		(1/4" ID-250' Roll) 9126-250

Custom hose assembly configurator is available at: <http://hoseconfig.precisionmedical.com>

		9130
Crimping Tool		
		
Hose Retractor		
1/4" Single Clamp	HR-1	
1/4" Double Clamp	HR-2	
5/16" Single Clamp	HR-3	
5/16" Double Clamp	HR-4	

Part Number	Autoclave	Method of Cleaning							
		Ethylene Oxide	Chlorine Dioxide	Hydrogen Peroxide	Gamma Ray	Electron Beam	Vaporized Hydrogen Peroxide (VHP)	Quaternary Cleaners	Sodium Hypochlorite "Bleach Solutions"
0154/0154C		X	X	X	X	X	X	X	X
0158/0158C		X	X	X	X	X	X	X	X
0159/0159C	X	X	X	X	X	X	X	X	X
0159K/0159KC	X	X	X	X	X	X	X	X	X
0115	X	X	X	X	X	X	X	X	X

* Please Remove O-Ring for Autoclaving *

* Clean all Tubing Nipples after each patient use *

		Oxygen	0154
Hand Tight Tubing Nipple		Oxygen with Chain	0154C
		Oxygen	0158
Wing Nut Tubing Nipple		Oxygen with Chain	0158C
		Oxygen	0159
Hex Nut Tubing Nipple		Oxygen with chain	0159C
		Oxygen	0159K
Knurled Nut Tubing Nipple		Oxygen with chain	0159KC
		Clear Yellow Black White Green	0119 0120 0121 0122 0124
Disposable Tubing Nipple		*Note: sold in packs of 25 and 50. Add -25 or -50 at end of part number for those quantities	
		Straight	9073
Tubing Nipple (for Vacuum Regulator)		90°	9090
		6" Chain with Closed Ring	8004
Chain Ring Assembly		6" Chain with Open Ring	8005

	1/8" NPT	1/4" NPT
 Pipe Tee	7030	7031
 Street Tee	7032	7033
 Female Tee	7034	7035
 Plug	7019	7017
 Drop Ear	Tee 7081 Elbow 7082	
 Elbow	7092	

	1/8" NPT	1/4" NPT
 Plug	7013	7012
 45° Street Elbow	7056	7057
 90° Elbow	7058	7059
 90° Street Elbow	7063	7065
 Master Flow Valve	7064	7045
 Cross	7036	7037

 Hex Nipple	1/8" NPT Male x 1/8" NPT Male 1/4" NPT Male x 1/4" NPT Male 1/8" NPT Male x 1/4" NPT Male	7020 7021 7022
 Connector	1/8" NPT Female x 1/8" NPT Male 1/8" NPT Female x 1/4" NPT Male 1/4" NPT Female x 1/8" NPT Male 1/4" NPT Female x 1/4" NPT Male 1/8" NPT Female x 1/8" NPT Male x 2.5" Long	1017 7024 7025 7026 8955
 Coupler	1/8" NPT Female x 1/8" NPT Female 1/4" NPT Female x 1/4" NPT Female 1/4" NPT Female x 1/8" NPT Female	7027 7028 7029
 Hose Barb Connector	1/8" NPT Male x 1/8" Diameter 1/8" NPT Male x 3/16" Diameter 1/8" NPT Male x 1/4" Diameter 1/4" NPT Male x 1/4" Diameter 1/8" NPT Male x 5/16" Diameter 1/8" NPT Female x 1/4" Diameter 1/4" NPT Female x 1/4" Diameter 1/4" NPT Female x 5/16" Diameter	7038 7039 7040 7043 8952 7046 7049 8951
 Swiveling Hose Barbs	1/8" NPT Male x 5/16" Diameter 1/8" NPT Male x 1/4" Diameter	507765 507764
 Pipe Nipple	1/8" NPT Male x 0.75" Long 1/8" NPT Male x 1.5" Long 1/8" NPT Male x 2" Long 1/8" NPT Male x 2.5" Long 1/4" NPT Male x 1.5" Long 1/4" NPT Male x 2" Long 1/4" NPT Male x 2.5" Long 1/8" NPT Male x 3" Long	7044 7050 7051 7052 7053 7054 7055 7062
 Cap	1/8" NPT Female	7018
 Hose Splicer	1/4" Inside Diameter	7041

Replacement Parts		
	Inlet DISS Female Oxygen x Outlet DISS Male Oxygen	PM15-4
	Inlet DISS Female Air x Outlet DISS Male Air	PM15-41
	Inlet DISS Female Air x Outlet DISS Male Air 90° Down	PM15-43
	Inlet DISS Female Oxygen 90° Up x Outlet DISS Male Oxygen 90° Down	PM15-44
	Inlet DISS Male Air 90° Down x Outlet DISS Female Air 90° Up	PM15-45
	Jar	507588
	Filter	507587
	Jar and Filter	1271

Neb Drives

		Oxygen	Air
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8126 8127 8128 8129	8226 8227 8228 8229
Ohmeda			
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8136 8137 8138 8139	8236 8237 8238 8239
Chemetron			
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8116 8117 8118 8119	8216 8217 8218 8219
DISS Hand Tight			
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8146 8147 8148 8149	8246 8247 8248 8249
Puritan-Bennett with Knob			
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8166 8167 8168 8169	8266 8267 8268 8269
Schrader			
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8156 8157 8158 8159	8256 8257 8258 8259
Oxequip with Knob			
	6 to 7 l/min 7 to 8 l/min 8 to 9 l/min 10 to 11 l/min	8176 8177 8178 8179	8276 8277 8278 8279
Medstar			

Rescue Flows

	Oxygen 14 to 16 l/min
	DISS Hand Tight 8115 Ohmeda 8125 Chemetron 8135 Puritan-Bennett 8145 Oxequip 8155 Schrader 8165 Medstar 8175
	Oxygen 24 to 26 l/min
	DISS Hand Tight 8114 Ohmeda 8124 Chemetron 8134 Puritan-Bennett 8144 Oxequip 8154 Schrader 8164 Medstar 8174

	Oxygen	Air	Vacuum	Nitrous Oxide
FRENCH (AFNOR)				
 Probe x 1/8" NPT Male	505979	505980	505981	
 Probe x 1/4" Hose Barb	505766	505767	505768	

EUROPEAN (NIST)				
 Nut/Nipple x 1/8" NPT Male	505752	505772	505773	
 Nut/Nipple x 1/4" Hose Barb	505774	505775	505776	
 Body Adaptor x 1/8" NPT Male	505777	505778	505779	
 Body Adaptor x 1/4" Hose Barb	505780	505781	505782	

AUSTRALIAN (SIS)				
 Nut/Nipple x 1/8" NPT Male	505769	505770	505771	

	Oxygen	Air	Vacuum	Nitrous Oxide
SCANDINAVIAN (AGA)				
 Probe x 1/8" NPT Male	505763	505764	505765	
 Probe x Hose Barb Swivel	1/4" Hose Barb 505763S	1/4" Hose Barb 505764S	1/4" Hose Barb 505765S	

GERMAN (DIN)				
 Probe x 1/8" NPT Male	505753	505761	505762	
 Probe x Hose Barb Swivel	1/4" Hose Barb 505753S	1/4" Hose Barb 505761S	1/4" Hose Barb 505762S	

BRITISH (BOC)				
 Probe x 1/8" NPT Male	6180	6184	6182	6183
 Probe x Hose Barb Swivel	1/4" Hose Barb 6180S	1/4" Hose Barb 6184S	1/4" Hose Barb 6182S	1/4" Hose Barb 6183S

 PK Schrader	6130	6230		
--	------	------	--	--



PrecisionMedical

The Difference is Precision.

300 Held Drive | Northampton, PA 18067

precisionmedical.com

T: 800-272-7285 | F: 800-352-1240