

TRUDELL
HEALTHCARE SOLUTIONS

vyaireTM
MEDICAL



Multi-LinkTM X2 ECG System

Disposable single-patient-use leadwires

The **Multi-LinkTM X2** single-patient use leadwire system is a fully disposable and versatile solution.

The leadwires remain with the patient throughout their patient journey, even when transferred between departments.

Disposable leadwires that move with the patient reduce cross-contamination and help lower the rate of hospital acquired infections (HAI).

Finding the right connection in your healthcare facility is no longer a concern as the Multi-LinkTM X2 ECG System is compatible with the following brands:

- GE[®] Healthcare
- Philips[®]
- Spacelabs[®]
- Nihon Kohden[®]
- Mindray[®]


Multi-Link X2 ECG System SKUs, AHA

Extra-long leadwire configuration				
Monitor input	SKU	Monitor brand	Leads	Cable length
	2090100-041	GE	3 or 5	9 in/23 cm
	2090100-043		6	
	2090100-045		12	
	2090100-051	Philips	3 or 5	9 in/23 cm
	2090100-053		6	
	2090100-055		12	
	2090100-077	Spacelabs	3	9 in/23 cm
	2090100-071		5	
	2090100-075		12	
	2090100-081	Nihon Kohden	3	9 in/23 cm
	2090100-083		6	
	2090100-061	Mindray	3 or 5	9 in/23 cm
	2090100-065		12	
SKU				
	SKU	Connection style	Leads	Lead length
	2052133-405	Grabber	3	12 ft/3.6 m
	2052133-406		5	
	2052133-408		V	
	2086367-006		6	
	2052104-405	Snap	3	12 ft/3.6 m
	2052104-406		5	
	2052104-408		V	
	2086349-006		6	

Traditional configuration				
Monitor input	SKU	Monitor brand	Leads	Cable length
	2017003-001	GE	3 or 5	12 ft/3.6 m
	2017005-001		6	
	2017006-001		12	
	2090100-052	Philips	3 or 5	12 ft/3.6 m
	2090100-054		6	
	2090100-056		12	
	2090100-078	Spacelabs	3	12 ft/3.6 m
	2090100-072		5	
	2090100-076		12	
	2090100-082	Nihon Kohden	3	12 ft/3.6 m
	2090100-084		6	
	2090100-062	Mindray	3 or 5	12 ft/3.6 m
	2090100-066		12	
SKU		Connection style	Leads	Lead length
	2052133-005	Grabber	3	39 in/100 cm
	2052133-006		5	
	2086367-003		6	
	2052133-008		V	
	2052133-027		3	
	2052133-007		5	
	2086367-004	Snap	6	39 in/100 cm
	2052104-005		3	
	2052104-006		5	
	2086349-003		6	
	2052104-008		V	
	2052104-027		3	
	2052104-007		5	29 in/75 cm
	2086349-004		6	

**For inquiries, please contact
your Territory Manager or:**

1-800-265-5495 customerservice@trudellhs.com



TRUDELL
HEALTHCARE SOLUTIONS