

Oxygen Monitor and AccuO2 Oxygen Analyzer



*"Maximize your blender capability with
Precision Medical's modular add-on accessories"*



The Precision Medical Oxygen Monitor and AccuO2 Oxygen Analyzer have been expertly engineered to meet the demands of critical settings. Using precise engineering principles, the PM5900 and PM5950 provide users gas analysis capabilities with easy-to-use, low-maintenance features. Both units are backed by Precision Medical's 2-year warranty.

AccuO2 Oxygen Analyzer PM5950



Lightweight

Compact and lightweight design.


Simple

Easy one touch calibration.

Long Battery Life

5,000 + hours operation on 2 AA batteries.

Oxygen Monitor PM5900



Intelligent

The only Oxygen monitor that retains settings when batteries are removed.

Unique

Three standard alarms including high and low, check sensor and low battery.

Durable

Smart designed high impact housing for unmatched performance under the most demanding conditions.

Versatile

Engineered with simple mounting brackets, soft side grips and nonskid support feet.

User Friendly

Easy-to-read screen and simple key pad operation.

Accurate

Maintains +/- 3 percent accuracy.

Ordering Info and Accessories



Oxygen Monitor

Designed to accurately view oxygen concentrations including settable high and low alarms

Order Information:

PM5900 - includes Oxygen Monitor with Coil Cable, Diverter, Sensor, T Piece and 4 AA batteries



AccuO2 Oxygen Analyzer

Designed for accurate, reliable and low maintenance measurement of volume percent oxygen

Order Information:

PM5950 - includes Oxygen Analyzer with Coil Cable, Diverter, Sensor and T Piece,



SensO2 Oxygen Sensor Supply

Specifically designed for the Precision Medical Line of blenders, attaches to the left Primary outlet

Order Information:

PM5300SO



High Flow Helper

Versatile triple tubing connection adapter

Order Information:

1MFAHF



Sensor Tee

Easily connects any oxygen sensor to any flowmeter

Order Information:

1MFATP



Monitor / Analyzer Mounting Bracket

Durable, lightweight bracket for greater usability

Order Information:

508420

Ordering Info and Accessories



Coil Cable

Replacement sensor cable for PM5900 and PM5950

Order Information:

504937



505013

505014

Pole Mount

Dovetail pole mount

Order Information:

505013 Vertical Pole Mount

505014 Horizontal Pole Mount



Monitor V Mount Bracket

Rear dovetail replacement for PM5900 and PM5950

Order Information:

505010 Monitor V Mount Bracket

505120 Screw - V Mount Bracket



Wall Mount

Dovetail wall mount

Order Information:

505189



T Piece

Single use Sensor Tee piece

Order Information:

505126



Oxygen Sensor

Replacement galvanic oxygen sensor for PM5900 and PM5950

Order Information:

504877



Diverter

Sensor diverter

Order Information:

505344