

PowerVac+ Aspirator



PM63

The PowerVac+ Aspirator is the perfect solution for all your suction needs. Designed for continuous and intermittent use, this powerful and durable suction pump combines the performance expected for an institution with the ease of operation needed for homecare use. It features a 0-200 mmHg adjustable vacuum regulator gauge and delivers maximum vacuum of up to 250 mmHg. PowerVac+ is a simple-to-use stationary aspirator with a durable design to withstand the rigors of hospital and homecare use. It features a hospital grade vacuum regulator, an anti-bacteria filter design to prevent contamination, uses oil-less, permanent lubrication, and includes a collection canister with overflow protection.

Features

- Continuous and intermittent operation
- Hospital grade two-mode vacuum regulator
- Oil-free, permanent lubrication
- Movable canister holder for placement on either side
- Collection canister with overflow protection
- Easy-to-clean and infection control-friendly ABS plastic cabinet
- Replaceable inline hydrophobic bacteria filter

Product Specifications

- Vacuum range:**
Adjustable: 0-200 mmHg
Fully adjusted: ≥ 250 mmHg
- Vacuum modes:**
Regulated, intermittent & off
- Intermittent cycle time:**
Factory set at approximately 16 seconds "on" and 8 seconds "off"
- Height:** 10.1 in (25.7 cm)
- Width:** 7.4 in (18.7 cm)
- Length:** 15.9 in (40.3 cm)
- Weight:** 21.2 lbs (9.6 kg)
- Suction canister:** 2000cc
(Disposable with float shut-off)

Warranty

- 3-year on device

Ordering Information

- PowerVac+ Aspirator (PM63)**
Includes PowerVac+ with 2000cc canister, 6' patient suction tubing, 10" suction canister tubing, inline hydrophobic bacteria filter, plastic elbow connector and hospital-grade power cord

Aspirator Accessories

Canisters

502519-10

800cc canister with lid (case of 10) - PM60 and PM65

1906-12

1200cc canister with lid (case of 12) - PM61 and PM63

502687-12

2000cc canister with lid (case of 12) - PM61 and PM63

502540

Canister bracket - PM61 and PM63

Roll Stands

PM60-1

Roll stand for use with any Precision Medical aspirator or compressor

PM60-2

Roll stand with canister bracket for use with any Precision Medical aspirator

Power Supply and Cord

503974

Battery charger power supply - PM65

502954

Auto DC power cord - PM65

Carry Bag

502956

Carry bag - PM65

Tubing and Filters

1690

Inline hydrophobic bacteria filter - PM60 and PM63

1955

Patient suction tubing (6 feet)

502012

Inline hydrophobic bacteria filter assembly - PM60

502391

Small air filter - PM63

502471

Large air filter - PM63

502690

Inline hydrophobic bacteria filter & tube assembly - PM61 and PM63

502698

Plastic tubing replacement kit (2" and a 12" section) - PM60

503293

Inline hydrophobic bacteria filter assembly with tubing and elbow - PM65

506047

Inline hydrophobic bacteria filter - PM65

