

PowerVac Aspirator

Features

Continuous operation

Automatic resetting thermal limit switch prevents overheating

Dual pump protection systems with float shutoff and inline hydrophobic filter safeguards

Easy-to-clean and infection control-friendly ABS plastic cabinet

Replaceable inline hydrophobic bacteria filter

Product Specifications

Vacuum range: 0-622 mmHg

Vacuum mode: Continuous

Height: 10.1 in (25.7 cm)

Width: 7.4 in (18.7 cm)

Length: 15.9 in (40.3 cm)

Weight: 20.2 lbs (9.2 kg)

Suction canister: 1200cc
(Disposable with float shut-off)

Warranty

5-year on device

Ordering Information

PowerVac Aspirator (PM61)

Includes PowerVac with 2000cc canister, 6' patient suction tube, 10" suction canister tubing, inline hydrophobic bacteria filter and hospital-grade power cord



PM61

PowerVac Aspirator is built around a high-performance suction pump that delivers up to 622 mmHg of vacuum. The wide range of suction allows for operation in homecare and hospital settings. PowerVac is designed for continuous use and has a built-in safety system that protects the device from overheating. It features a dual protection system that includes a float shut-off in the canister lid and an inline hydrophobic bacteria filter to prevent back flow of fluids and airborne contaminants from entering the pump. PowerVac also has a unique, movable canister holder which allows the caregiver to rotate the canister to either side of the unit for better placement near the point of care.

Aspirator Accessories

Canisters

502519-10

800cc canister with lid (case of 10) - PM60 and PM65

1906-12

1200cc canister with lid (case of 12) - PM61 and PM63

502687-12

2000cc canister with lid (case of 12) - PM61 and PM63

502540

Canister bracket - PM61 and PM63

Roll Stands

PM60-1

Roll stand for use with any Precision Medical aspirator or compressor

PM60-2

Roll stand with canister bracket for use with any Precision Medical aspirator

Power Supply and Cord

503974

Battery charger power supply - PM65

502954

Auto DC power cord - PM65

Carry Bag

502956

Carry bag - PM65

Tubing and Filters

1690

Inline hydrophobic bacteria filter - PM60 and PM63

1955

Patient suction tubing (6 feet)

502012

Inline hydrophobic bacteria filter assembly - PM60

502391

Small air filter - PM63

502471

Large air filter - PM63

502690

Inline hydrophobic bacteria filter & tube assembly - PM61 and PM63

502698

Plastic tubing replacement kit (2" and a 12" section) - PM60

503293

Inline hydrophobic bacteria filter assembly with tubing and elbow - PM65

506047

Inline hydrophobic bacteria filter - PM65

