

SafeSet™

Blood Sampling System


Reduce blood waste and
increase the efficiency of
the blood sampling process.



A closed, needlefree in-line
blood sampling system that
reduces blood waste, while
minimizing the risk of IV line
contamination and transmission
of infectious disease.

icumedical
human connections

In response to the problems and risks associated with blood management, ICU Medical has developed the SafeSet Blood Sampling System. The SafeSet is an intuitive, easy-to-use system that helps conserve blood while protecting against exposure to bloodborne pathogens and increasing the efficiency and safety of the blood sampling process.



Safety and Efficiency of the SafeSet System Defined

- › SafeSet has a closed, in-line design that protects patient's blood from exposure to bloodstream infections and protects the caregiver by limiting exposure to bloodborne pathogens.
- › Reduces blood waste because the discard or clearing volume can be reinfused, which may reduce the risk of iatrogenic anemia.
- › Two needlefree access options protect the clinician from exposure to infectious diseases.
- › Impenetrable polycarbonate case protects sheath and plunger from exposure to non-sterile environment.
- › Kit configurations to address both adult and the neonatal/pediatric patient population.
- › Available as a standalone system, or with the ICU Medical Transpac® IV Pressure Monitoring Kit, or as an add-on to other manufacturers pressure monitoring systems.

PRODUCT SPECIFICATIONS

Pressure Tubing	Internal Volume, .055 mL/in
Reservoir	Adult, 12 mL graduated to 10 mL Neonatal/pediatric 3 mL graduated
Without Latex Components	Yes

PRODUCT ACCESSORIES

Shielded Blunt Cannula	42303-02
Shielded Blunt Cannula with Tube Holder	42327-02

SafeSet, a blood conservation system to help you reach your blood management goals.

icumedical
human connections

800.824.7890 | www.icumed.com