



**Up to 97%
savings in waste
management
costs**

Serres Nemo

**SERRES SUCTION BAG EMPTYING DEVICE FOR SAFER
AND MORE SUSTAINABLE WASTE MANAGEMENT**

**Cost
effective**

**ECO-
friendly**


Safe

Convenient

serres

serres

Cut costs and CO₂ emissions with Serres Nemo



SERRES NEMO is a device for the disposal of surgical fluid waste collected in suction bags. It helps hospitals cut costs and CO₂ emissions. Serres Nemo provides a safer working environment for the operating room and waste management personnel.

Cost effective

- Up to 97% savings in waste management costs

Eco-friendly

- Up to 97% savings in CO₂ emissions due to reduced amount of waste

Safe

- Improved occupational safety
- Minimised risk of contamination


Convenient

- Minimal amount of waste to manoeuvre
- Device fits into a tight space



“We gained 80% cost savings in hazardous fluid waste management.”

Heiko Sahm
Waste Manager and
Occupational Safety Specialist
Klinikum Gütersloh, Germany



An easy way to empty Serres Suction Bags

What makes Serres Nemo an eco-friendly solution? - Significantly lower CO₂ emissions

Most suction bags are incinerated when full. But using Serres Nemo to empty the bags before transportation and incineration, helps cut down emissions significantly.

TRANSPORTATION

100 km



12g CO₂

Full 2 litre suction bags



1g CO₂

Emptied with Serres Nemo

WASTE INCINERATION



150g CO₂


Full 2 litre suction bags



5g CO₂

Emptied with Serres Nemo

97%
less CO₂



“We don’t need to move around with full and heavy suction bags anymore; Serres Nemo is accessible to everyone.”

Veronique Fabry
OR Nurse, Clinique Beau Soleil
Montpellier, France

Convenient, safe, and easy to use:

1. Insert the full Serres Suction Bag into the sealed container.
2. Close the lid with a safe click.
3. Press the start button.
4. Dispose of the empty Serres Suction Bag in accordance with hospital waste disposal guidelines.



Compact device: 330 x 515 x 1060 mm

 Easy to clean with Serres washing tray.



The leader in sustainable surgical fluid waste management

Serres is a progressive leader in developing sustainable surgical fluid management solutions.

We support over 60,000 operations a day, in more than 50 markets, ensuring that healthcare professionals can focus on patient care.

Responsible product design and co-creation with our customers have resulted in unmatched high-quality products with only 1-in-a-million product failures.

Our passion is to provide eco-friendly solutions that enable our customers to reach their sustainability targets.



We report based on GRI standards
Find our sustainability report from [Serres.com](https://www.serres.com)

serres

