

TRU**TEST**® Arterial Blood Gas Sampling Kits

TRU**TEST**® offers a full range of arterial blood gas sampling kits. They are used to draw arterial blood samples for in vitro diagnostics, pH blood gas and selected whole blood electrolyte testing.

- Vented and non-vented syringes
- Dry balanced lithium-zinc heparin
- Thin-walled needle technology
- 3cc and 1cc volumes
- Sterile
- Latex free

All kits are compatible with a large range of blood gas analyzers used in clinical settings.



Vented Syringe Options

Item Number: T222532320

UOM: Case/100

Components:

- 3cc vented syringe, luer lock
- 25u lithium-zinc heparin
- 23G needle
- Safe Tip safety mechanism and cap



Item Number: T22151000

UOM: Case/100

Components:

- 1cc vented syringe, luer slip
- 15u lithium-zinc heparin
- Syringe cap (*not shown below*)



Item Number: T222532320M

UOM: Case/100

Components:

- 25u lithium-zinc heparin
- 3cc vented syringe, luer lock
- 23G needle
- Safe Tip safety mechanism and cap
- Alcohol and iodine prep pads
- Gauze and bandaid
- Biohazard ice bag
- Patient label



Cord Gas Kit

Item Number: T22251020

UOM: Case/100

Components:

- Two 1cc non-vented syringes, luer lock
- 25u lithium/zinc heparin
- Two 22G needles
- Safe Tip safety mechanisms and cap
- Biohazard ice bag
- Venous and arterial labels



Non-Vented Syringe Options

Item Number: T22253000

UOM: Case/100

Components:

- 3cc non-vented syringe, luer lock
- 25u lithium-zinc heparin
- Syringe cap (*not shown below*)



Item Number: T22251000

UOM: Case/100

Components:

- 1cc non-vented syringe, luer slip
- 25u lithium-zinc heparin
- Syringe cap (*not shown below*)



Item Number: T22251001

UOM: Case/100

Components:

- 1cc non-vented syringe, luer slip
- 3cc non-vented waste syringe, luer lock
- 25u lithium-zinc heparin
- Syringe caps (*not shown below*)



Item Number: T222535601

UOM: Case/100

Components:

- 3cc non-vented syringe, luer lock
- 5-6cc non-vented waste syringe, luer lock
- 25u lithium-zinc heparin
- Syringe caps (*not shown below*)



Convenience Packs

TRUtest® offers a variety of ABG convenience pack configurations that will streamline your clinical practice.

- Save time (*no product assembly*)
- Reduce costs (*all products ordered under one code*)
- Reduce packaging waste

M023400MU	3cc Syringe, 25u, 23g, Vented Safe Tip w/CHG-Alc Wipe, Gauze, Ice Bag
M024023EC10106	3cc Syringe, 25u, 23g, Vented Safe Tip w/CHG-Alc Wipe, Gauze, Bandaid, Patient Label
M024023RHEC	3cc Syringe, 25u, 23g, Vented Safe Tip w/Alc Pad, Gauze, Bandaid, Bio Bag
M02P4023LTRUBPMU	3cc Syringe, 25u, 23g, Vented Safe Tip w/Alc Pad, Gauze, Bandaid, Patient Label, Ice Bag
M02P4023LTRUMU	3cc Syringe, 25u, 23g, Vented Safe Tip w/Alc Pad, Gauze, Ice Bag
M02P9000RHST	1cc Syringe, 25u Non-Vented w/Red Dot Sticker

For more information, please contact your local THS Representative.
customers@trudellhs.com | 1-800-265-5495