

## Flexible face masks

Our high-quality disposable, flexible face masks offer a soft face plate, teardrop shape and superior grip. Pediatric sizes are offered in an optional child-friendly bubble gum scent. Easy to inflate or deflate air cushion with top valve style. Vyaire offers color-coded packaging for easy size identification.



Cat. no.	Description	Qty.
<b>Flex dome, no valve, air filled cushion, clear</b>		
6800EU	Premie size 0	20
<b>Flex dome, top valve, air filled cushion, clear</b>		
6810EU	Neonate size 1	20
6820	Infant size 2	20
6821	Infant size 2	100
6830	Toddler size 3	20
6831	Toddler size 3	100
6840	Child size 4	20
6841	Child size 4	100
6850	Adult size 5	20
6851	Adult size 5	100
6860	XL adult size 6	20
6861	XL adult size 6	100
<b>Flex dome, top valve, air filled cushion, clear</b>		
6820BG	Infant size 2, bubble gum scented	20
6821BG	Infant size 2, bubble gum scented	100
6830BG	Toddler size 3, bubble gum scented	20
6831BG	Toddler size 3, bubble gum scented	100
6840BG	Child size 4, bubble gum scented	20
6841BG	Child size 4, bubble gum scented	100

## Round face masks

Our traditional line of non-flexible, rounded dome face masks provides several styles to suit individual needs. These disposable masks feature a dual inflation valve (tail valve) that permits cushion adjustment with a syringe or by mouth. Soft, pliable, air-filled cushion molds easily to face to deliver a low-pressure seal. Vyair masks are made without latex.



Cat. no.	Description	Qty.
<b>Round rigid dome, air filled cushion, clear</b>		
5210	Neonate size 1, w/o valve	20
5220	Infant size 2, w/o valve	20
5221EU	Infant size 2, w/o valve	100
5230	Toddler size 3, tail valve	20
5231EU	Toddler size 3, tail valve	100
5230X	Toddler size 3, w/o valve	20
5240	Child size 4, tail valve	20
5241	Child size 4, tail valve	100
5250EU	Adult size 5, tail valve	20
5251VS	Adult size 5, tail valve	100
5211	Adult size 5, tail valve	100
<b>Rigid dome, air filled cushion, clear</b>		
5246	Adult size 5, tail valve	100
5245	Adult size 5, tail valve	20

## Standard face masks

See the difference in our newest line of standard disposable face masks - featuring color-coded hook rings, brilliant clarity and semi-soft teardrop face plates.



Cat. no.	Description	Qty.
<b>Standard, top valve, air filled cushion, clear</b>		
4810VS	Neonate size 1, red ring	20
4820VS	Infant size 2, green ring	20
4830VS	Toddler size 3, purple ring	20
4840VS	Child size 4, yellow ring	20
4852VS	Small adult size 5, white ring	20
4860VS	Medium adult size 6, blue ring	20
4861VS	Medium adult size 6, blue ring	100
4870VS	Large adult size 7, orange ring	20