



## Container Repair and Maintenance


Trudell Surgical Instrument Care Services is committed to preserving the life of your surgical instruments and containers used to store and transport these surgical assets.

Container integrity plays a vital role in ensuring that the sterility of your instruments is not compromised. Containers and wire baskets should be included in preventative maintenance programs where they will be inspected, repaired and cleaned as required.

### **Repair and Refurbishment of Aesculap STERILCONTAINER™ Systems:**

- Preventative maintenance extends the useful life of a container by repairing or replacing damaged parts and professionally inspecting the container to assure optimal performance.
- Refurbishing containers at the same time the set is being repaired reduces set downtime.
- TSICS will remove materials such as ink, tape and marker from container surfaces without damaging the outer anodized surface. This value added process is provided free of charge as part of any container repair service.





## The Complete Solution

Our commitment to the repair and refurbishment of surgical instruments is demonstrated by the recent re-design of a 2500 ft<sup>2</sup> area in the head office dedicated exclusively to TSCIS and surgical instrument repair. TSCIS adheres to strict process and performance standards established by the manufacturer. By honouring these stringent processes customers are guaranteed the highest quality workmanship in the industry.

Our technicians are able to repair and refurbish all makes of surgical instruments. Regardless of brand, only new factory replacement parts are used to maintain the integrity of the instrument.

Canadian hospitals spend a large portion of their operating budget on surgical instruments. Our commitment to quality workmanship can help protect the longevity of that investment by:

- Reducing replacement costs.
- Minimizing set down time by scheduling routine repairs.
- Maintaining the highest performance level of all instruments via traditional craftsmanship.



Trudell Surgical Instrument Care Services